



Cytokeratin, LMW

USA IVD

This is used to distinguish poorly differentiated carcinomas from non-epithelial neoplasms.

Mouse Monoclonal, Clone: AE1

Control: Prostate. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁ **Cellular localization:** Cytoplasmic



Prostate

Catalog no.	0.1 ml concentrated
YL-2010-1	0.5 ml concentrated
YL-2010-2	1ml concentrated
YL-2010-3	1ml Prediluted
YL-2010-4	7 ml Prediluted
YL-2010-5	5 Positive control slides
YL-2010	