




Desmin

USA IVD

This is used to demonstrate the myogenic components and derivation of tumors.

Mouse Monoclonal, Clone: D33

Control: Muscle. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Cytoplasmic



Muscle

Catalog no.	Volume
De-2010-1	0.1 ml concentrated
De-2010-2	0.5 ml concentrated
De-2010-3	1ml concentrated
De-2010-4	1ml Prediluted
De-2010-5	7 ml Prediluted
De-2010	5 Positive control slides

