




## E-Cadherin

USA IVD

This is useful in distinguishing adenocarcinoma from mesothelioma.

**Rabbit Monoclonal, Clone: EP700Y**

**Control:** Breast. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG **Cellular localization:** Membranous



Breast

<b>Catalog no.</b>	<b>Volume</b>
EC-2010-1	0.1 ml concentrated
EC-2010-2	0.5 ml concentrated
EC-2010-3	1ml concentrated
EC-2010-4	1ml Prediluted
EC-2010-5	7 ml Prediluted
EC-2010	5 Positive control slides

