



# Epithelial Specific Antigen/Ep-CAM (MOC-31)

USA IVD

This is used to distinguish adenocarcinoma from mesothelioma and hepatocellular carcinoma.

**Mouse Monoclonal, Clone: MOC-31**

**Control:** Colon Carcinoma. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub>/K **Cellular localization:** Cytoplasmic



Colon Carcinoma

## Catalog no.

	<b>Volume</b>
ES-2010-1	
ES-2010-2	0.1 ml concentrated
ES-2010-3	0.5 ml concentrated
ES-2010-4	1ml concentrated
ES-2010-5	1ml Prediluted
ES-2010	7 ml Prediluted
	5 Positive control slides

