




# Epstein Barr Virus

USA ASR

This targets the 60kDa latent membrane protein encoded by the BNLF1 gene of the virus.

**Mouse Monoclonal, Clone: CS1-4**

**Control:** EBV Infected Tissue. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub> **Cellular localization:** Cyto plasmic



EBV Infected Tissue

Catalog no.	Volume
EB-2010-1	0.1 ml concentrated
EB-2010-2	0.5 ml concentrated
EB-2010-3	1ml concentrated
EB-2010-4	1ml Prediluted
EB-2010-5	7 ml Prediluted
EB-2010	5 Positive control slides