




FLI-1

USA IV D

This is used in the pathogenesis of ES/pPNET

Mouse Monoclonal, Clone: MRQ-1

Control: PNET. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG_{2b} **Cellular localization:** Nuclear



PNET

Catalog no.	Volume
FL-2010-1	0.1 ml concentrated
FL-2010-2	0.5 ml concentrated
FL-2010-3	1ml concentrated
FL-2010-4	1ml Prediluted
FL-2010-5	7 ml Prediluted
FL-2010	5 Positive control slides