




## FSH

USA IVD

This is a marker in classification of pituitary tumors and the study of pituitary disease.

### Rabbit Polyclonal

**Control:** Pituitary. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Cellular localization:** Cytoplasmic



Pituitary

<b>Catalog no.</b>	0.1 ml concentrated
FS-2010-1	0.5 ml concentrated
FS-2010-2	1ml concentrated
FS-2010-3	1ml Prediluted
FS-2010-4	7 ml Prediluted
FS-2010-5	5 Positive control slides
FS-2010	

Page 96

### Volume