




Factor VIII

USA IV D

This helped establish the endothelial nature of some lesions of disputed histogenesis.

Rabbit Polyclonal

Control: Placenta. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Cytoplasmic



Placenta

Catalog no.	0.1 ml concentrated
FA-2010-1	0.5 ml concentrated
FA-2010-2	1ml concentrated
FA-2010-3	1ml Prediluted
FA-2010-4	7 ml Prediluted
FA-2010-5	5 Positive control slides

Page 90

Volume