



# Fascin


USA IVD

This is a prognostic marker in neuroendocrine neoplasms of the lung and ovarian cancer.

**Mouse Monoclonal, Clone: 55k-2**

**Control:** Hodgkin's Lymphoma. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C.

**Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub> **Cellular localization:** Cytoplasmic



Hodgkin's Lymphoma

| <b>Catalog no.</b> | <b>Volume</b>             |
|--------------------|---------------------------|
| Fc-2010-1          | 0.1 ml concentrated       |
| Fc-2010-2          | 0.5 ml concentrated       |
| Fc-2010-3          | 1ml concentrated          |
| Fc-2010-4          | 1ml Prediluted            |
| Fc-2010-5          | 7 ml Prediluted           |
| Fc-2010            | 5 Positive control slides |