



Fascin


USA IVD

This is a prognostic marker in neuroendocrine neoplasms of the lung and ovarian cancer.

Mouse Monoclonal, Clone: 55k-2

Control: Hodgkin's Lymphoma. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C.

Shelf life: Up to 36 months. **Isotype:** IgG₁ **Cellular localization:** Cytoplasmic



Hodgkin's Lymphoma

Catalog no.	Volume
Fc-2010-1	0.1 ml concentrated
Fc-2010-2	0.5 ml concentrated
Fc-2010-3	1ml concentrated
Fc-2010-4	1ml Prediluted
Fc-2010-5	7 ml Prediluted
Fc-2010	5 Positive control slides