




Follicular Dendritic Cell

USA IVD

This is useful in the identification of follicular dendritic cell matrix found in normal lymph nodes and tonsil tissue.

Mouse Monoclonal, Clone: CNA.42

Control: Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgM/K **Cellular localization:** Cytoplasmic



Tonsil

Catalog no.	Volume
FD-2010-1	0.1 ml concentrated
FD-2010-2	0.5 ml concentrated
FD-2010-3	1ml concentrated
FD-2010-4	1ml Prediluted
FD-2010-5	7 ml Prediluted
FD-2010	5 Positive control slides