




# Follicular Dendritic Cell

USA IVD

This is useful in the identification of follicular dendritic cell matrix found in normal lymph nodes and tonsil tissue.

**Mouse Monoclonal, Clone: CNA.42**

**Control:** Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgM/K **Cellular localization:** Cytoplasmic



Tonsil

<b>Catalog no.</b>	<b>Volume</b>
FD-2010-1	0.1 ml concentrated
FD-2010-2	0.5 ml concentrated
FD-2010-3	1ml concentrated
FD-2010-4	1ml Prediluted
FD-2010-5	7 ml Prediluted
FD-2010	5 Positive control slides