




Glial Fibrillary Acidic Protein (GFAP)

USA IVD

This is used to distinguish neoplasms of astrocytic origin from other neoplasms in the central nervous system.

Rabbit Monoclonal, Clone: EP672Y

Control: Brain **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** EP672Y, IgG, G-A-5, IgG₁ **Cellular localization:** Cytoplasmic



Brain

Catalog no.

GF-2010-1

GF-2010-2

GF-2010-3

GF-2010-4

GF-2010-5

GF-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides