




# Hepatitis B Virus Surface Antigen

USA ASR

This labels the cell cytoplasm infected with hepatitis B virus.

**Mouse Monoclonal, Clone: S1-210**

**Control:** Hepatitis B Infected Tissue **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C.  
**Shelf life:** Up to 36 months. **Isotype:** IgG<sub>2b</sub> **Cellular localization:** Cytoplasmic.



Hepatitis B Infected Tissue

<b>Catalog no.</b>	<b>Volume</b>
HE-2010-1	0.1 ml concentrated
HE-2010-2	0.5 ml concentrated
HE-2010-3	1ml concentrated
HE-2010-4	1ml Prediluted
HE-2010-5	7 ml Prediluted
HE-2010-6	5 Positive control slides.