




Hepatitis B Virus Core Antigen

USA ASR

This labels cell nuclei infected with hepatitis B virus.

Rabbit Polyclonal

Control: Hepatitis B Infected Tissue **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C.
Shelf life: Up to 36 months. **Cellular localization:** Nuclear.



Hepatitis B Infected Tissue

Catalog no.

HB-2010-1

HB-2010-2

HB-2010-3

HB-2010-4

HB-2010-5

HB-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides.