




# Hepatocyte Specific Antigen

USA IVD

This is useful in differentiating hepatocellular carcinomas with adenoid features from adenocarcinomas.

**Mouse Monoclonal, Clone: OCH1E5**

**Control:** Liver **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub>/K **Cellular localization:** Cytoplasmic



Liver

## Catalog no.

HS-2010-1

HS-2010-2

HS-2010-3

HS-2010-4

HS-2010-5

HS-2010

## Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides.