




Herpes Simplex Virus I & II

USA ASR

Rabbit Polyclonal

Control: HSV Infected Tissue **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Cellular localization:** Cytoplasmic, Nuclear



HSV Infected Tissue

Catalog no.

HV-2010-1

HV-2010-2

HV-2010-3

HV-2010-4

HV-2010-5

HV-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides.