




# Human Chorionic Gonadotropin (hCG)

USA IVD

This antibody detects; cells and tumors of trophoblastic origin such as choriocarcinoma.

## Rabbit Polyclonal

**Control:** Placenta **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. | **Cellular localization:** Cytoplasmic



Placenta

## Catalog no.

HC-2010-1

HC-2010-2

HC-2010-3

HC-2010-4

HC-2010-5

HC-2010

5 Positive control slides

## Volume

0.1 ml, concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted