




INI-1

USA IVD

This is a marker of medulloblastomas and sPNEts.

Mouse Monoclonal, Clone: MRQ-27

Control: Brain **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Isotype:** IgG_{2a} **Cellular localization:** Nuclear



Brain

Catalog no.

IN-2010-1

IN-2010-2

IN-2010-3

IN-2010-4

IN-2010-5

IN-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides.