




# IgD

USA IVD

This is useful when identifying leukemias, plasmacytomas, and B-cell lineage derived lymphomas.

## Rabbit Polyclonal

**Control:** Tonsil **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Cellular localization:** Cytoplasmic



Tonsil

## Catalog no.

ID-2010-1

ID-2010-2

ID-2010-3

ID-2010-4

ID-2010-5

ID-2010

5 Positive control slides .

## Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted