




Kappa

USA IVD

This detects surface immunoglobulin on normal and neoplastic B-cells.

Mouse Monoclonal, Clone: L1C1

Control: Tonsil **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Cytoplasmic



Tonsil

Catalog no.

KA-2010-1

KA-2010-2

KA-2010-3

KA-2010-4

KA-2010-5

KA-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides.