




Ki-67

USA IVD

This is a good marker of neoplasms such as grade II astrocytoma, oligodendroglioma, colon carcinoma, and breast carcinoma.

Rabbit Monoclonal, Clone: SP6

Control: Tonsil **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Nuclear



Tonsil

Catalog no.

KI-2010-1

KI-2010-2

KI-2010-3

KI-2010-4

KI-2010-5

KI-2010

7 ml Prediluted

5 Positive control slides.

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted