



LH

USA IVD

This is a marker in classification of pituitary tumors and the study of pituitary disease.

Rabbit Polyclonal

Control: Pituitary **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Cellular localization:** Cytoplasmic



Pituitary

Catalog no.

LH-2010-1

LH-2010-2

LH-2010-3

LH-2010-4

LH-2010-5

LH-2010

7 ml Prediluted

5 Positive control slides.

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted