




# Lambda

USA IVD

This detects surface immunoglobulin on normal and neoplastic B-cells.

**Mouse Monoclonal, Clone: Lamb14**

**Control:** Tonsil **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>2a</sub> **Cellular localization:** Cytoplasmic



Tonsil

## Catalog no.

LA-2010-1

LA-2010-2

LA-2010-3

LA-2010-4

LA-2010-5

LA-2010

## Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides.