




# Lysozyme

USA IVD

This may identify histiocytic neoplasia, large lymphocytes, and classifying lymphoproliferative disorders.

## Rabbit Polyclonal

**Control:** Tonsil **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Cellular localization:** Cytoplasmic



Tonsil

## Catalog no.

LY-2010-1

LY-2010-2

LY-2010-3

LY-2010-4

LY-2010-5

LY-2010

1ml Prediluted

7 ml Prediluted

5 Positive control slides.

## Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated