




MSH6

USA IVD

Mutations in this gene have been identified in individuals with hereditary nonpolyposis colon cancer. In *E. coli*, the MutS protein helps in the recognition of mismatched nucleotides, prior to their repair.

Mouse Monoclonal, Clone: 44

Control: Colon Embedded: Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG1 **Cellular localization:** Nuclear



Colon

Catalog no.

MH-2010-1

MH-2010-2

MH-2010-3

MH-2010-4

MH-2010-5

MH-2010

7 ml Prediluted

5 Positive control slides.

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted