




# MUC1

USA IVD

This is in apical surfaces of most epithelial cells in breast, GI, respiratory, and GU tracts.

**Mouse Monoclonal, Clone: MRQ-17**

**Control:** Breast **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub> **Cellular localization:** Cytoplasmic



Breast

## Catalog no.

MU-2010-1

MU-2010-2

MU-2010-3

MU-2010-4

MU-2010-5

MU-2010

## Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides.