




# MUC2

USA IVD

This is specifically expressed in goblet cells of the small intestine and colon.

**Mouse Monoclonal, Clone: MRQ-18**

**Control:** Colon **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub>/K **Cellular localization:** Cytoplasmic



Colon

## Catalog no.

MU2-2010-1

MU2-2010-2

MU2-2010-3

MU2-2010-4

MU2-2010-5

MU2-2010

## Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides .