




MUC5AC

USA IVD

This is preferentially expressed in the normal stomach and respiratory tract; esophageal carcinomas-67%, gastric carcinomas-58%.

Mouse Monoclonal, Clone: MRQ-19

Control: Stomach **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Isotype:** IgG1 **Cellular localization:** Cytoplasmic



Stomach

Catalog no.

MU5-2010-1

MU5-2010-2

MU5-2010-3

MU5-2010-4

MU5-2010-5

MU5-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides .