




MUC6

USA IVD

Mouse Monoclonal, Clone: MRQ-20

Control: Stomach **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁ **Cellular localization:** Cytoplasmic



Stomach

Catalog no.

MU6-2010-1

MU6-2010-2

MU6-2010-3

MU6-2010-4

MU6-2010-5

MU6-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides.