




Melanoma

USA IVD

This stains fetal and neonatal melanocytes, junctional and blue nevus cells, and malignant melanocytes.

Mouse Monoclonal, Clone: HMB-45

Control: Melanoma **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Cytoplasmic



Melanoma

Catalog no.

ME-2010-1

ME-2010-2

ME-2010-3

ME-2010-4

ME-2010-5

ME-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides.