




Myelin Basic Protein

USA IVD

The immunoreactivity for this in granular cell tumors strengthens the concept of a Schwann cell derivation of lesions.

Rabbit Polyclonal

Control: Brain **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Cellular localization:** Cytoplasmic



Brain

Catalog no.

MP-2010-1

MP-2010-2

MP-2010-3

MP-2010-4

MP-2010-5

MP-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides .