




# Myogenin

USA IVD

This labels the nuclei of myoblasts in developing muscle tissue.

**Mouse Monoclonal, Clone: F5D**

**Control:** Rhabdomyosarcoma **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub>/K **Cellular localization:** Nuclear



Rhabdomyosarcoma

## Catalog no.

MN-2010-1

MN-2010-2

MN-2010-3

MN-2010-4

MN-2010-5

MN-2010

## Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides.