



Myoglobin

USA IVD

This provides a specific, sensitive, and practical procedure for the identification of tumors of muscle origin.

Rabbit Polyclonal

Control: Skeletal Muscle **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG1/K **Cellular localization:** Cytoplasmic



Skeletal Muscle

Catalog no.

MB-210-1

MB-2010-2

MB-2010-3

MB-2010-4

MB-2010-5

MB-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides