




# Myosin, Smooth Muscle

USA IVD

This helps to distinguish between benign sclerosing breast lesions and carcinomas.

**Mouse Monoclonal, Clone: SMMS-1**

**Control:** Breast **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub>/K **Cellular localization:** Cytoplasmic



Breast

## Catalog no.

MO-210-1

MO-2010-2

MO-2010-3

MO-2010-4

MO-2010-5

MO-2010

## Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides