




## NSE

USA IVD

Along with other antibodies, this can help detect neural and neuroendocrine lineage.

**Mouse Monoclonal, Clone: E27**

**Control:** Pancreas **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub> **Cellular localization:** Cytoplasmic



Pancreas

### Catalog no.

NS-2010-1

NS-2010-2

NS-2010-3

NS-2010-4

NS-2010-5

NS-2010

### Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides