




Napsin A

USA IVD

This helps to distinguish primary lung adenocarcinomas from adenocarcinomas of other organs.

Rabbit Polyclonal

Control: Lung Adenocarcinoma, Kidney **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Cellular localization:** Cytoplasmic



Lung Adenocarcinoma

Catalog no.

NA-2010-1

NA-2010-2

NA-2010-3

NA-2010-4

NA-2010-5

NA-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides .