




Neurofilament

USA IVD

This stains an antigen localized in a number of neural, neuroendocrine, and endocrine tumors.

Mouse Monoclonal, 2F11

Control: Brain **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Cytoplasmic



Brain

Catalog no.

NF-210-1

NF-2010-2

NF-2010-3

NF-2010-4

NF-2010-5

NF-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides.