




## Oct-2

### USA IVD

This is a transcription factor of the POU homeo domain family that regulates B-cell specific genes.

#### Mouse Monoclonal, Clone: MRQ-2

**Control:** Tonsil **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub>/K **Cellular localization:** Nuclear



Tonsil

#### Catalog no.

OC-2010-1

OC-2010-2

OC-2010-3

OC-2010-4

OC-2010-5

OC-2010

#### Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides .