




Oct-4

USA IVD

This is a transcription factor that maintains and regulates pluripotency in embryonic stem and germ cells.

Mouse Monoclonal, Clone: MRQ-10

Control: Seminoma **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Isotype:** IgG **Cellular localization:** Nuclear



Seminoma

Catalog no.

OT-2010-1

OT-2010-2

OT-2010-3

OT-2010-4

OT-2010-5

OT-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides