




PAX-5

USA IVD

This is a marker for B-cells including B-lymphoblastic neoplasms and maturation stage.

Mouse Monoclonal, Clone: 24

Control: Tonsil **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** 24, IgG₁, SP34, IgG **Cellular localization:** Nuclear



Tonsil

Catalog no.

PX-2010-1

PX-2010-2

PX-2010-3

PX-2010-4

PX-2010-5

PX-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides