




# PD-1

USA IVD

This is a marker for angioimmunoblastic lymphoma.

**Mouse Monoclonal, Clone: MRQ-22**

**Control:** Tonsil **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub> **Cellular localization:** Cytoplasmic



Tonsil

## Catalog no.

PD-2010-1

PD-2010-2

PD-2010-3

PD-2010-4

PD-2010-5

PD-2010-6

## Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides