




# PLAP

USA IV D

This can discriminate between germ cell tumors and other neoplasms.

**Mouse Monoclonal, Clone: NB10**

**Control:** Placenta **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** NB10, IgG/K SP15, IgG. **Cellular localization:** Nuclear



Placenta

## Catalog no.

PL-2010-1

PL-2010-2

PL-2010-3

PL-2010-4

PL-2010-5

PL-2010

## Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides