




# PMS2

USA IVD

**Mouse Monoclonal, Clone: MRQ-28**

**Control:** Colon **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8<sup>0</sup>C. **Shelf life:** Up to 36 months. **Cellular localization:** Nuclear



Colon

## Catalog no.

PM-2010-1

PM-2010-2

PM-2010-3

PM-2010-4

PM-2010-5

PM-2010

## Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml prediluted

5 Positive control slides