




# PSAP

USA IVD

This is helpful in finding the site of origin in cases of metastatic carcinoma of the prostate.

**Mouse Monoclonal, Clone: PASE/4LJ**

**Control:** Prostate **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub> **Cellular localization:** Cytoplasmic



Prostate

## Catalog no.

PP-2010-1

PP-2010-2

PP-2010-3

PP-2010-4

PP-2010-5

PP-2010

## Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides