



# Parvovirus B19


USA ASR

This targets the capsid proteins VP1 and VP2 on human parvovirus.

**Mouse Monoclonal, Clone: R92F6**

**Control:** Parvovirus Infected Tissue **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C.

**Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub> **Cellular localization:** Cytoplasmic, Nuclear



Parvovirus Infected Tissue

## Catalog no.

PV-2010-1

PV-2010-2

PV-2010-3

PV-2010-4

PV-2010-5

PV-2010

## Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides .