



Pneumocystis jiroveci (carinii)


USA ASR

This reacts with an epitope on an organism which is resistant to formalin, picric acid, paraffin, alcohol, and xylene.

Mouse Monoclonal, Clone: 3F6

Control: P. Jiroveci Infected Tissue **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C.

Shelf life: Up to 36 months. **Isotype:** IgM/K **Cellular localization:** Membranous



P. Jiroveci Infected Tissue

Catalog no.

PJ-2010-1

PJ-2010-2

PJ-2010-3

PJ-2010-4

PJ-2010-5

PJ-2010

Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides