




Prolactin

USA IVD

This is a marker of pituitary tumors and the study of pituitary disease.

Rabbit Polyclonal

Control: Pituitary **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Isotype:** **Cellular localization:** Cytoplasmic



Pituitary

Catalog no.

PT-2010-1

PT-2010-2

PT-2010-3

PT-2010-4

PT-2010-5

PT-2010

Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides