




# Spectrin

USA IVD

This is useful in the diagnosis of erythroid leukemias.

**Mouse Monoclonal, Clone: RBC2/3D5**

**Control:** Bone Marrow **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>2b</sub>/K **Cellular localization:** Membranous



Bone Marrow

## Catalog no.

SN-2010-1

SN-2010-2

SN-2010-3

SN-2010-4

SN-2010-5

SN-2010

## Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides