




Synaptophysin

USA IV D

This reacts with neuroendocrine cells of human adrenal medulla, carotid body, skin, pituitary, thyroid, lung, pancreas, and gastrointestinal mucosa.

Rabbit Polyclonal

Control: Pancreas **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Cellular localization:** Cytoplasmic



Pancreas

Catalog no.

SY-2010-1

SY-2010-2

SY-2010-3

SY-2010-4

SY-2010-5

SY-2010

Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides